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f_{up} : fraction unbound in plasma. R_b : blood to plasma ratio. CL: clearance. CL_{int} : intrinsic clearance. CL_h : hepatic clearance. F_{renal} : the fraction of drug excreted by renal elimination. a) Chan *et al.*⁴⁾ b) Stringer *et al.*¹⁵⁾ c) Dave *et al.*¹⁶⁾ d) Obach *et al.*¹⁷⁾ e) Uchimura *et al.*¹⁸⁾ f) R_b was assumed as 1. g) These values were obtained by dividing the reference value (% Urinary Excretion) by one hundred. * Ratio = in vitro Predicted CL_h /in vivo Observed CL.

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1. S.M. Paul, D.S. Mytelka, A.L. Schacht, et al. How to improve R&D productivity: the pharmaceutical industry's grand challenge. *Nat. Rev. Drug Discov.* 2010, 9, 203-214.
2. I. Kola, J. Landis. Can the pharmaceutical industry reduce attrition rates? *Nat. Rev. Drug Discov.* 2004, 3, 711-715.
3. J. Arrowsmith, P. Miller. Phase II and Phase III attrition rates 2011-2012. *Nat. Rev. Drug Discov.* 2013, 12, 569.
4. C. Rodriguez-Antona, M.T. Donato, M.J. Gómez-Lechón, et al. Cytochrome P450 expression in human hepatocytes and hepatoma cell lines: molecular mechanisms that determine lower expression in cultured cells. *Xenobiotica.* 2002, 32, 6, 505-520.
5. L. Di, P. Trapa, L. M. Tremaine, et al. A Novel Relay Method for Determining Low-Clearance Values. *Drug Metab. Dispos.* 2012, 40, 9, 1860-1865.
6. S.R. Khetani, S.N. Bhatia. Microscale culture of human liver cells for drug development. *Nat. Biotechnol.* 2008, 26, 1, 120-126.
7. T.S. Chan, H. Yu, D. Tweedie, et al. Meeting the Challenge of Predicting Hepatic Clearance of Compounds Slowly Metabolized by Cytochrome P450 Using a Novel Hepatocyte Model. *HepatoPac. Drug Metab. Dispos.* 2019, 47, 1, 58-66.
8. I. Hultman, C. Vedin, M. Darnell, et al. Use of HμREL® human co-culture system for prediction of intrinsic clearance and metabolite formation for slowly metabolized compounds. *Mol. Pharmaceutics.* 2016, 13, 8, 2796-2807.
9. T. Ohkura, K. Ohta, T. Ogihara, et al. Evaluation of Human Hepatocytes Cultured by Three-dimensional Spheroid Systems for Drug Metabolism. *Drug Metab. Pharmacokinet.* 2014, 29, 5, 373-378.
10. C.C. Bell, D.F.G. Hendriks, M. Ingelman-Sundberg, et al. Characterization of primary human hepatocyte spheroids as a model system for drug-induced liver injury, liver function and disease. *Sci. Rep.* 2016, 6, 25187.
11. S.U. Vorrink, Y. Zhou, V.M. Lauschke, et al. Prediction of Drug-Induced Hepatotoxicity Using Long-Term Stable Primary Hepatic 3D Spheroid Cultures in Chemically Defined Conditions. *Toxicol. Sci.* 2018, 163, 2, 655-665.
12. T. Ogihara, H. Arakawa, H. Kojima, et al. Utility of human hepatocyte spheroids without feeder cells for evaluation of hepatotoxicity. *J. Toxicol. Sci.* 2017, 42, 4, 499-507.
13. T. Takebe, K. Sekine, T. Ogaeri, et al. Ran-Ran Zhang, Yasuharu Ueno, Yun-Wen Zheng, Naoto Koike, Shinsuke Aoyama, Yasuhisa Adachi, and Hideki Taniguchi. Vascularized and functional human liver from an iPSC-derived organ bud transplant. *Nature.* 2013, 499, 481-484.
14. T. Shinozawa, H. Yoshikawa, T. Takebe. Reverse Engineering Liver Buds Through Self-Driven Condensation and Organization Towards Medical Application. *Developmental Biology.* 2016, 420, 221-229.
15. A.R. Baudy, M.A. Otieno, W.R. Proctor, et al. Liver microphysiological systems development guidelines for safety risk assessment in the pharmaceutical industry. *Lab Chip.* 2020, 20, 215-225.
16. C. Yamasaki, C. Tateno, K. Yoshizato, et al. Growth and differentiation of colony-forming human hepatocytes in vitro. *J. Hepatol.* 2006, 44, 749-757.
17. T. Takezawa, K. Ozaki, T. Shimo-Oka, et al. Collagen Vitrigel: A Novel Scaffold That Can Facilitate a Three-Dimensional Culture for Reconstructing Organoids. *Cell Transplant.* 2004, 13, 463-473.
18. A. Oshikata-Miyazaki, T. Takezawa. Development of an oxygenation culture method for activating the liver-specific functions of HepG2 cells utilizing a collagen vitrigel membrane chamber. *Cytotechnology.* 2016, 68, 5, 1801-1811.
19. S. Sanoh, Y. Naritomi, S. Ohta, et al. Predictability of plasma concentration-time curves in humans using single-species allometric scaling of chimeric mice with humanized liver. *Xenobiotica.* 2015, 45, 7, 605-614.
20. M. Miyamoto, S. Iwasaki, H. Hirabayashi, et al. Comparison of predictability for human pharmacokinetics parameters among monkeys, rats, and chimeric mice with humanised liver. *Xenobiotica.* 2017, 47, 12, 1052-1063.
21. C. Tateno, K. Yoshizato. Long-term cultivation of adult rat hepatocytes that undergo multiple cell divisions and express normal parenchymal phenotypes. *Am J Pathol* 1996;148:383-392.
22. R. Watari, M. Kakiki, K. Kusano, et al. A long-term culture system based on a collagen vitrigel membrane chamber that supports liver-specific functions of hepatocytes isolated from mice with humanized livers. *J. Toxicol. Sci.* 2018, 43, 8, 521-529.
23. R. Watari, M. Kakiki, K. Kusano, et al. Prediction of Human Hepatic Clearance for Cytochrome P450 Substrates via a New Culture Method Using the Collagen Vitrigel Membrane Chamber and Fresh Hepatocytes Isolated from Liver Humanized Mice. *Biol. Pharm. Bull.* 2019, 42, 3, 348-353.
24. F. Da-Silva, X. Boulenc, P. Poulin, et al. Improving Prediction of Metabolic Clearance Using Quantitative Extrapolation of Results Obtained From Human Hepatic Micropatterned Cocultures Model and by Considering the Impact of Albumin Binding. *J. Pharm. Sci.* 2018, 107, 1957-1972.