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1. 竹村 晃典, 伊藤 晃成. In vitro 肝毒性評価の進歩. 日本薬理学雑誌. 2020, 155, 6, 401-405.
2. K. Harada, H. Kohara, T. Yukawa, K. Matsumiya, and T. Shinozawa. Cell-based high-throughput screening for the evaluation of reactive metabolite formation potential. *Toxicol In Vitro*. 2021, 74, 105159.
3. Y. Sanuki, T. Araki, O. Nakazono, and K. Tsurui. A rapid mitochondrial toxicity assay utilizing rapidly changing cell energy metabolism. *J Toxicol Sci*. 2017, 42, 3, 349-358.
4. H. Kohara, P. Bajaj, K. Yamanaka, A. Miyawaki, K. Harada, K. Miyamoto, T. Matsui, Y. Okai, M. Wagoner, and T. Shinozawa. High-Throughput Screening to Evaluate Inhibition of Bile Acid Transporters Using Human Hepatocytes Isolated From Chimeric Mice. *Toxicol Sci*. 2020, 173, 2, 347-361.
5. T. Koga, K. Takeuchi, and K. Umehara. Establishment of an in vitro cholestasis risk assessment system using two-dimensional cultured HepaRG cells and 12 bile acids. *J Toxicol Sci*. 2023, 48, 1, 47-56.
6. M. Kawaguchi, T. Nukaga, S. Sekine, A. Takemura, T. Susukida, S. Oeda, A. Kodama, M. Hirota, H. Kouzuki, and K. Ito. Mechanism-based integrated assay systems for the prediction of drug-induced liver injury. *Toxicol Appl Pharmacol*. 2020, 394, 114958.
7. Y. Ikeyama, T. Sato, A. Takemura, S. Sekine, and K. Ito. Hypoxia/reoxygenation exacerbates drug-induced cytotoxicity by opening mitochondrial permeability transition pore: Possible application for toxicity screening. *Toxicol In Vitro*. 2020, 67, 104889.
8. A. Takemura, S. Gong, T. Sato, M. Kawaguchi, S. Sekine, Y. Kazuki, T. Horie, and K. Ito. Evaluation of Parent- and Metabolite-Induced Mitochondrial Toxicities Using CYP-Introduced HepG2 cells. *J Pharm Sci*. 2021, 110, 9, 3306-3312.
9. A. Takemura, S. Ishii, Y. Ikeyama, K. Esashika, J. Takahashi, and K. Ito. New in vitro screening system to detect drug-induced liver injury using a culture plate with low drug sorption and high oxygen permeability. *Drug Metab Pharmacokinet*. 2023, 52, 100511.
10. J.A. Dykens and Y. Will. The significance of mitochondrial toxicity testing in drug development. *Drug Discov Today*. 2007, 12, 17-18, 777-785.
11. C. Liu, S. Sekine, and K. Ito. Assessment of mitochondrial dysfunction-related, drug-induced hepatotoxicity in primary rat hepatocytes. *Toxicol Appl Pharmacol*. 2016, 302, 23-30.
12. G. Labbe, D. Pessayre, and B. Fromenty. Drug-induced liver injury through mitochondrial dysfunction: mechanisms and detection during preclinical safety studies. *Fundam Clin Pharmacol*. 2008, 22, 4, 335-353.
13. Y. Ikeyama, T. Sato, A. Takemura, S. Sekine, and K. Ito. Successful energy shift from glycolysis to mitochondrial oxidative phosphorylation in freshly isolated hepatocytes from humanized mice liver. *Toxicol In Vitro*. 2020, 65, 104785.
14. Y. Sakai, M. Nishikawa, F. Evenou, M. Hamon, H. Huang, K.P. Montagne, N. Kojima, T. Fujii, and T. Niino. Engineering of implantable liver tissues. *Methods Mol Biol*. 2012, 826, 189-216.
15. C. Funk, C. Ponelle, G. Scheuermann, and M. Pantze. Cholestatic potential of troglitazone as a possible factor contributing to troglitazone-induced hepatotoxicity: in vivo and in vitro interaction at the canalicular bile salt export pump (Bsep) in the rat. *Mol Pharmacol*. 2001, 59, 3, 627-635.
16. T. Susukida, S. Sekine, M. Nozaki, M. Tokizono, and K. Ito. Prediction of the Clinical Risk of Drug-Induced Cholestatic Liver Injury Using an In Vitro Sandwich Cultured Hepatocyte Assay. *Drug Metab Dispos*. 2015, 43, 11, 1760-1768.
17. T. Shinozawa, M. Kimura, Y. Cai, N. Saiki, Y. Yoneyama, R. Ouchi, H. Koike, M. Maezawa, R. R. Zhang, A. Dunn, A. Ferguson, S. Togo, K. Lewis, W. L. Thompson, A. Asai, and T. Takebe. High-Fidelity Drug-Induced Liver Injury Screen Using Human Pluripotent Stem Cell-Derived Organoids. *Gastroenterology*. 2021, 160, 3, 831-846.
18. S. U. Vorrink, Y. Zhou, M. Ingelman-Sundberg, and V. M. Lauschke. Prediction of Drug-Induced Hepatotoxicity Using Long-Term Stable Primary Hepatic 3D Spheroid Cultures in Chemically Defined Conditions. *Toxicol Sci*. 2018, 163, 2, 655-665.
19. C. C. Bell, A. C. A. Dankers, V. M. Lauschke, R. Sison-Young, R. Jenkins, C. Rowe, C.E. Goldring, K. Park, S. L. Regan, T. Walker, C. Schofield, A. Baze, A. J. Foster, D. P. Williams, A. W. M. van de Ven, F. Jacobs, J. V. Houdt, T. Lähteenmäki, J. Snoeys, S. Juhila, L. Richert, and M. Ingelman-Sundberg. Comparison of Hepatic 2D Sandwich Cultures and 3D Spheroids for Long-term Toxicity Applications: A Multicenter Study. *Toxicol Sci*. 2018, 162, 2, 655-666.
20. 横井 毅. 薬物性肝障害評価研究の進歩. 日本薬理学雑誌. 2020, 155, 5, 323-328.
21. J. Uetrecht. Mechanistic Studies of Idiosyncratic DILI: Clinical Implications. *Front Pharmacol*. 2019, 10, 837.
22. R. Kato and J. Uetrecht. Supernatant from Hepatocyte Cultures with Drugs That Cause Idiosyncratic Liver Injury Activates Macrophage Inflammasomes. *Chem Res Toxicol*. 2017, 30, 6, 1327-1332.
23. R. Kato, Y. Ijiri, T. Hayashi, and J. Uetrecht. The 2-Hydroxyiminostilbene Metabolite of Carbamazepine or the Supernatant from Incubation of Hepatocytes with Carbamazepine Activates Inflammasomes: Implications for Carbamazepine-Induced Hypersensitivity Reactions. *Drug Metab Dispos*. 2019, 47, 10, 1093-1096.
24. W. R. Proctor, A. J. Foster, J. Vogt, C. Summers, B. Middleton, M. A. Pilling, D. Shienson, M. Kijanska, S. Ströbel, J. M. Kelm, P. Morgan, S. Messner, and D. Williams. Utility of spherical human liver microtissues for prediction of clinical drug-induced liver injury. *Arch Toxicol*. 2017, 91, 8, 2849-2863.
25. C. C. Bell, B. Chouhan, L. C. Andersson, H. Andersson, J. W. Dear, D. P. Williams, and M. Söderberg. Functionality of primary hepatic non-parenchymal cells in a 3D spheroid model and contribution to acetaminophen hepatotoxicity. *Arch Toxicol*. 2020, 94, 4, 1251-1263.
26. A. Segovia-Zafra, M. Villanueva-Paz, A. S. Serras, G. Matilla-Cabello, A. Bodoque-García, D. E. Di Zeo-Sánchez, H. Niu, I. Álvarez-Álvarez, L. Sanz-Villanueva, S. Godec, I. Milisav, P. Bagnaninchi, R. J. Andrade, M. I. Lucena, J. C. Fernández-Checa, F. J. Cubero, J. P. Miranda, and L. J. Nelson. Control compounds for preclinical drug-induced liver injury assessment: Consensus-driven systematic review by the ProEuroDILI network. *J Hepatol*. 2024, 81, 4, 630-640.
27. M. Arai, O. Yokosuka, K. Fujiwara, H. Kojima, T. Kanda, H. Hirasawa, and H. Saisho. Fulminant hepatic failure associated with benzbromarone treatment: a case report. *J Gastroenterol Hepatol*. 2002, 17, 5, 625-626.
28. A. Takemura, S. Gong, S. Sekine, and K. Ito. Inhibition of biliary network reconstruction by benzbromarone delays recovery from pre-existing liver injury. *Toxicology*. 2019, 423, 32-41.
29. W. Yang, Y. Tao, Y. Wu, X. Zhao, W. Ye, D. Zhao, L. Fu, C. Tian, J. Yang, F. He, and L. Tang. Neutrophils promote the development of reparative macrophages mediated by ROS to orchestrate liver repair. *Nat Commun*. 2019, 10, 1, 1076.
30. N. Tanimizu, N. Ichinohe, H. Suzuki, and T. Mitaka. Prolonged oxidative stress and delayed tissue repair exacerbate acetaminophen-induced liver injury in aged mice. *Aging (Albany NY)*. 2020, 12, 19, 18907-18927.

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