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Table 1. Water Absorption Capacities of Various Absorbents

Sample Description		Water absorption Capacity (g) per gram of dried sample
Test Group	SP ^{a)}	1349
	FP ^{a)}	1295.4
Control Group	Pulp	3.8
	Silica Gel	1.4
	Ion-Exchange Resin	2.5
	PVA	4.6
	High-Grade Synthetic High-polymer Absorbent ^{b)}	249.4
	Anionic Synthetic High-polymer Absorbent ^{c)}	363.6

a) Bioabsorbent samples were produced using different culture conditions, i.e., by changing the culture medium's carbon sources; SP: Sucrose, FP: Fructose

b) High Grade Synthetic High polymer Absorbent : poly acrylate / PYA derivative (co-polymer of acrylate and vinyl alcohol.)

c) Anionic Synthetic High-polymer Absorbent : Polyacrylamide derivative, MW. 350 X 10⁴

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
Fig. 1. Moisture retention capacities of various absorbents. Conditions: 33% relative humidity, 20°C, in desiccator containing a saturated of magnesium chloride.

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
Fig. 2. Bioabsorbent and Synthetic High Polymer Absorbent Water Absorption Capacities in the Presence of Various NaCl Concentrations. ○, bioabsorbent (SP sample); ●, high-grade water absorbing synthetic high-polymer (Sumika Gel S-50).

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
Fig. 4. HPLC Charts for Standard Samples and the Hydrolyzate of the Bioabsorbent.

A: Standard samples. 1. Glucuronic acid 2. Glucose 3. Mannose 4. Rhamnose 5. Fucose.
 B: Hydrolyzate of polysaccharide component including high and low-molecular-weight biopolymer before separation.
 C: Hydrolyzate of the high-molecular-weight acid biopolymer (the bioabsorbent).
 D: Hydrolyzate of the low-molecular-weight neutral biopolymer.

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Table 2. Effects of amino acid and yeast extract on biopolymer production (viscosity : cps)

Amino acid	Viscosity (cps)	Amino acid	Viscosity (cps)
Yeast extract	2,140	Valine	900
Methionine	1,200	Tyrosine	920
Leucine	120	Histidine	540
Isoleucine	200	Proline	20
Phenylalanine	540	Aspartic acid	40
Tryptophan	32	Glycine	3,400
Threonine	480	Citrulline	4,000
Glutamic acid	2,656	Lysine	520
Alanine	320	Arginine	2,440
Serine	240	Cysteine	320
Ornithine	800	none	—

— : No viscosity was found,

Rotary viscometer, no.3 rotor, 30rpm.

Culture conditions : Glucose 1%, urea 0.1%, vitamin mixture, metal mixture, pH7.2, Erlenmeyer flask (300ml), rotary shaker, 180rpm/min, 7d at 30°C.

Vitamin mixture : p-aminobenzoic acid 20 µg, inositol 0.2mg, vitamin A 2 µg, nicotinic acid 40 µg, vitamin B₁₂ 0.5µg, vitamin B₁ 40 µg, Ca-pantothenate 0.4 mg, vitamin B₂ 40 µg, vitamin B₆ 40µg / 1,000ml of the medium.

Metal mixture : FeSO₄·7H₂O 5mg, Na₂WO₄·2H₂O 0.5 mg, Na₂MoO₄·2H₂O 0.5 mg, MnSO₄·5H₂O 0.5 mg, CaCl₂ 10mg / 1,000 ml of the medium.

Table 3. Viscosity comparison of various thickening agents

	Apparent viscosity (mPa·s)		
	30°C	50°C	78°C
B-16 polymer (untreated)	50	45	45
Xanthan gum	56	36	11
Wellan gum	46	40	20
Rhamzan gum	76	70	65
B-16 polymer (pH 12.25, 5 min)	175	168	147
B-16 polymer (pH 12.25, 80 min)	515	380	230

The viscosities were measured with a BM-type rotational viscometer (no.2 rotior, 30rpm, 1 min) after storing each polymer aqueous solution (0.1wt%) overnight at 30°C, 50°C and 78°C.

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<The viscosity vs. temperature>


Fig. 5. Apparent Viscosities vs. Temperature of the bioabsorbent from *A.latus* B-16 and other polymers.

●; B-16 polymer (0.2%), △; Xanthan gum (0.5%), □; Sodium alginate (1.0%), ◇; HPMC (1.0%).

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