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1. EURLMB, EURLMB SOP for the analysis of Paralytic shellfish toxins (PST) by precolumn HPLC-FLD according to OMA AOAC 2005.06, https://www.aesan.gob.es/en/CRLMB/docs/docs/metodos_analiticos_de_desarrollo/PST_OMA-AOAC2005-06_HPLC-FLD_EURLMB.pdf (2025年3月5日アクセス)
2. A. Turner, R. Hatfield, M. Rapkova-Dhanji, D. Norton, M. Algoet, D. Lees, Single-Laboratory Validation of a Refined AOAC HPLC Method 2005.06 for Oysters, Cockles, and Clams in UK Shellfish, *J. AOAC int.*, 2010, 93, pp.1482-1493.
3. Y. Oshima, Post-column derivatization HPLC method for the analysis of Paralytic shellfish toxins, *J. AOAC int.*, 1995, 78, pp528-532.
4. E. Hignutt, Jr., Suitability of Postcolumn Oxidation Liquid Chromatography Method AOAC 2011.02 for Monitoring Paralytic Shellfish Toxins in Alaskan Shellfish—Initial Pilot Study versus Mouse Bioassay and In-House Validation. *J. AOAC int.*, 2014, 97, pp.293-298.
5. M. J. Boundy, A. I. Selwood, D. T. Harwood, P. S. McNabb, A. D. Turner, Development of a sensitive and selective liquid chromatography-mass spectrometry method for high throughput analysis of paralytic shellfish toxins using graphitised carbon solid phase extraction, *J. Chromatogr. A*, 2015, 1384, p1-12
6. R. Watanabe, M. Kanamori, H. Yoshida, Y. Okumura, H. Uchida, R. Matsushima, H. Oikawa, T. Suzuki, Development of Ultra-Performance Liquid Chromatography with Post-Column Fluorescent Derivatization for the Rapid Detection of Saxitoxin Analogues and Analysis of Bivalve Monitoring Samples. *Toxins*, 2019, 11, p573.
7. (a) R. Watanabe, T. Suzuki, Y. Oshima, Preparation of Calibration Standards of N1-H Paralytic Shellfish Toxin Analogues by Large-Scale Culture of Cyanobacterium *Anabaena circinalis* (TA04). *Marine Drugs*, 2011, 9(3), 466-477. (b) T. Suzuki, R. Watanabe, A. Yoshino, H. Oikawa, H. Uchida, R. Matsushima, S. Nagai, T. Kamiyama, T. Yamazaki, M. Kawaguchi, T. Yaritha, A. Takatsu. Preparation of Diarrhetic Shellfish Toxins (DSTs) and Paralytic Shellfish Toxins (PSTs) by Large Algal Culture and Chemical Conversion. A. Lincoln MacKenzie [Ed]. *Marine and Freshwater Harmful Algae 2014. Proceedings of the 16th International Conference on Harmful*
8. O. Iwamoto, K. Nagasawa, Total Synthesis of (+)-Decarbamoylsaxitoxin and (+)-Gonyautoxin 3, *Org. Lett.* 2010, 12, pp2150-2153
9. N. Ohyabu, T. Nishikawa, M. Isobe. First asymmetric total synthesis of tetrodotoxin. *J. Am. Chem. Soc.* 2003. 125. pp.8798-8805
10. I. W. Burton, M. A. Quilliam, J. A. Walter, Quantitative ¹H NMR with External Standards: Use in Preparation of Calibration Solutions for Algal Toxins and Other Natural Products, *Anal. Chem.* 2005, 77 (10), 3123-3131
11. R. Watanabe, T. Suzuki, Y. Oshima, Development of quantitative NMR method with internal standard for the standard solutions of paralytic shellfish toxins and characterisation of gonyautoxin-5 and gonyautoxin-6, *Toxicon*, 2010, 56, pp. 589-595
12. R. Watanabe, C. Sugai, T. Yamazaki, R. Matsushima, H. Uchida, M. Matsumiya, A. Takatsu, T. Suzuki, Quantitative Nuclear Magnetic Resonance Spectroscopy Based on PULCON Methodology: Application to Quantification of Invaluable Marine Toxin, Okadaic Acid. *Toxins*, 2016, 8(10), 294.
13. 「qNMR プライマリーガイド」ワーキンググループ, qNMR プライマリーガイド 基礎から実践まで, 共立出版, 2015, p58
14. C. Y. Hong, Y. Kishi, Enantioselective total synthesis of (-)-decarbamoylsaxitoxin, *J. Am. Chem. Soc.* 1992, 114, pp. 7001-7006.
15. R. Watanabe, Y. Takayanagi, O. Chiba, S. Itoda, H. Ishizuka, M. Odagi, M. Ozawa, H. Uchida, R. Matsushima, K. Konoki, M. Yotsu-Yamashita, K. Nagasawa, T. Suzuki. Nontoxic Enantiomeric Reference Materials for Saxitoxins. *Anal. Chem.* 2022, 94, pp.11144-11150.
16. 経済産業省, 「化学兵器の禁止及び特定物質の規制等に関する法律」に係る特定物質の取扱いについて, https://www.meti.go.jp/policy/chemical_management/cwc/files/kaieihou_tokuteibushitsu_toriatsukai.pdf (2025年3月5日)
17. M. Kodama, S. Sato, T. Ogata, *Alexandrium tamarense* as a source of tetrodotoxin in the scallop *Patinopecten yessoensis*. Toxic phytoplankton blooms in the sea. 1993.
18. C. Dell'Aversano, J. A. Walter, I. W. Burton, D. J. Stirling, E. Fattorusso, M. A. Quilliam. Isolation and Structure Elucidation of New and Unusual Saxitoxin Analogues from Mussels. *J. Nat. Prod.* 2008, 71, pp. 1518-1523.
19. S. Numano, Y. Kudo, Y. Cho, K. Konoki, Y. Kaga, K. Nagasawa, M. Yotsu-Yamashita. Two new skeletal analogues of saxitoxin found in the scallop, *Patinopecten yessoensis*, as possible metabolites of paralytic shellfish toxins. *Chemosphere*, 2021, 278, p.130224.

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