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- Asada K, Inaba Y, Tateda-Suzuki E, Kuwahara-Arai K, Ito T, Hiramatsu K. Evolution and resistance expression of MRSA. Evaluation of beta-lactam antibiotics against a set of isogenic strains with different types of phenotypic expression. *Acta Biochemica Polonica*, **42** (4), 517-24 (1995).
- CDC. Staphylococcus aureus resistant to vancomycin --- United States, *Morbidity and Mortality Weekly Report*, **51**, 565-567 (2002).
- CDC. Vancomycin-resistant Staphylococcus aureus --- Pennsylvania, *Morbidity and Mortality Weekly Report*, **51**, 902 (2002).
- CDC. Vancomycin-Resistant Staphylococcus aureus --- New York, *Morbidity and Mortality Weekly Report*, **53**, 565-567 (2002).
- National Committee for Clinical Laboratory Standards. 2003. Methods for dilution antimicrobial susceptibility tests for bacteria that grow aerobically, 6th ed Approved standard NCCLS publication no. M7-A6. National Committee for Clinical Laboratory Standards, Wayne, Pa. (2003).
- Clinical and Laboratory Standards Institute, Methods for dilution antimicrobial susceptibility tests for bacteria that grow aerobically, 7th ed. Approved standard. CLSI publication no. M7- A7. Clinical and Laboratory Standards Institute, Wayne, Pa (2006)
- Fridkin SK, Hageman J, McDougal LK, Mohammed J, Jarvis WR, Perl TM, Tenover FC: Vancomycin-Intermediate Staphylococcus aureus Epidemiology Study Group. Epidemiological and Microbiological Characterization of Infections Caused by Staphylococcus aureus with Reduced Susceptibility to Vancomycin, United States, 1997-2001. *Clin Infect Dis.*, **15**, 36:429-39 (2003).
- Hiramatsu K, Aritaka N, Hanaki H, Kawasaki S, Hosoda Y, Hori S, Fukuchi Y, and Kobayashi I, Dissemination in Japanese hospitals of strains of Staphylococcus aureus of heterogeneously resistant to vancomycin, *Lancet*, **350**, 1668-71 (1997).
- Howe RA, Walsh TR. hGISA: seek and ye shall find. *Lancet*, **364**, 500-501 (2004).
- Ike Y, Arakawa Y, Ma X, Tatewaki K, Nagasawa M, Tomita H, Tanimoto K, Fujimoto S. Nationwide survey shows that methicillin-resistant Staphylococcus aureus strains heterogeneously and intermediately resistant to vancomycin are not disseminated throughout Japanese hospitals, *J Clin Microbiol.*, **39**:4445-51 (2001).
- Arakawa Y, Ike Y, Nagasawa M, Where has vancomycin-heterogeneously resistant *Staphylococcus aureus* gone? *Lancet*. **363**:1401 (2004).
- Liu H, Buescher G, Lewis N, Snyder S, Jungkind D, Detection of borderline oxacillin-resistant *Staphylococcus aureus* and differentiation from methicillin-resistant strains. *Eur J Clin Microbiol Infect Dis.*, **9** (10), 717-24 (1990).
- Sa-Leao R, Santos Sanches I, Dias D, Peres I, Barros RM, de Lencastre H. Detection of an archaic clone of Staphylococcus aureus with low-level resistance to methicillin in a pediatric hospital in Portugal and in international samples: relics of a formerly widely disseminated strain? *J Clin Microbiol.*, **37** (6), 1913-1920 (1999).
- Hanaki H, Yamaguchi Y, Yanagisawa C, Uehara K, Matsui H, Yamaguchi Y, Hososaka Y, Barada K, Sakai F, Itabashi Y, Ikeda S, Atsuda K, Tanaka H, Inamatsu T, Nagayama A, Sunakawa K. Investigation of beta-lactam antibiotic-induced vancomycin-resistant MRSA (BIVR). *J Infect Chemother.*, **11** (2), 104-106 (2005).
- Hanaki H, Yamaguchi Y, Nomura S, Haraga I, Nagayama A, Sunakawa K. Method of detecting beta-lactam antibiotic induced vancomycin resistant MRSA (BIVR). *Int J Antimicrob Agents.*, **23** (1), 1-5 (2004).
- Fey PD, Said-Salim B, Rupp ME, Hinrichs SH, Boxrud DJ, Davis CC, Kreiswirth BN, Schlievert PM. Comparative molecular analysis of community- or hospital-acquired methicillin-resistant *Staphylococcus aureus*, *Antimicrob Agents Chemother.*, **47** (1), 196-203 (2003).
- Storch GA, Rajagopalan L. Methicillin-resistant Staphylococcus aureus bacteremia in children, *Pediatr Infect Dis.*, **5** (1), 59-67 (1986).
- Jansen WT, Beitsma MM, Koeman CJ, van Wamel WJ, Verhoef J, Fluit AC. Novel mobile variants of staphylococcal cassette chromosome mec in *Staphylococcus aureus*. *Antimicrob Agents Chemother.*, **50** (6), 2072-2078 (2006).
- Takizawa Y, Taneike I, Nakagawa S, Oishi T, Nitahara Y, Iwakura N, Ozaki K, Takano M, Nakayama T, Yamamoto T. A Pantone-Valentine leucocidin (PVL) -positive community-acquired methicillin-resistant Staphylococcus aureus (MRSA) strain, another such strain carrying a multiple-drug resistance plasmid and other more-typical PVL-negative MRSA strains found in Japan. *J Clin. Microbiol.* **43** (7), 3356-3363 (2005).
- Peacock JE Jr, Moorman DR, Wenzel RP, Mandell GL. Methicillin-resistant Staphylococcus aureus: microbiologic characteristics, antimicrobial susceptibilities, and assessment of virulence of an epidemic strain, *J Infect Dis.*, **144** (6), 575-82 (1981).
- Hanaki H, Inaba Y, Sasaki K, Hiramatsu K. A novel method of detecting *Staphylococcus aureus* heterogeneously resistant to vancomycin (hetero-VRSA). *Jpn J Antibiot.*, **51** (8), 521-530 (1998).
- Hatano K, Yokota Y, Hanaki H, Sunakawa K. Combined effect of vancomycin or teicoplanin plus a beta-lactam antibiotic in mouse infection models caused by beta-lactam antibiotic-induced vancomycin resistant MRSA (BIVR) *Kansenshogaku Zasshi.*, **80** (3), 243-50 (2006).
- Haraga I, Nomura S, Fukamachi S, Ohjimi H, Hanaki H, Hiramatsu K, Nagayama A. Emergence of vancomycin resistance during therapy against methicillin-resistant Staphylococcus aureus in a burn patient--importance of low-level resistance to vancomycin. *Int J Infect Dis.*, **6** (4), 302-308 (2002).
- Takayama Y, Hanaki H, Irinoda K, Kokubun H, Yoshida K, Sunakawa K, Investigation of methicillin-resistant Staphylococcus aureus showing reduced vancomycin susceptibility isolated from a patient with infective endocarditis. *Int J Antimicrob Agents.*, **22** (6), 567-73 (2003).