

# Diagnostics Reagents List



- For clinical chemistry test and immunoassay
- Ready-to-use stable reagents
- Applicable to various clinical analyzers

# Clinical Chemistry

Test Item	Product Name	Method	Type	Package	Code	
<b>Enzyme</b>						
AST	CicaLiquid AST	JSCC	7170	Reagent 1	60 mL × 4	78213
				Reagent 2	20 mL × 4	77466
ALT	CicaLiquid ALT	JSCC	7170	Reagent 1	60 mL × 4	78215
				Reagent 2	20 mL × 4	77468
CK	CicaFit CK II	JSCC	70/20	Reagent 1	60 mL × 2	78551
				Reagent 2	18 mL × 2	
CK-MB	CicaLiquid CK-MB	Immunoinhibition	7170	Reagent 1	35 mL × 2	77822
				Reagent 2	10 mL × 2	77823
γ-GT	CicaLiquid γ-GT J	JSCC	7170	Reagent 1	60 mL × 4	77530
				Reagent 2	20 mL × 4	77531
LD	CicaFit LD-IFCC	IFCC	7170	Reagent 1	60 mL × 4	78325
				Reagent 2	33 mL × 4	78326
ALP	CicaLiquid ALP-IFCC	IFCC	7170	Reagent 1	55 mL × 4	78313
				Reagent 2	20 mL × 4	78314
AMY	CicaFit AMY-G7	JSCC (Et-G7-pNP)	7170	Reagent 1	60 mL × 4	78037
				Reagent 2	20 mL × 4	78038
p-AMY	CicaFit p-AMY-G7	Immunoinhibition (Et-G7-pNP)	SPS	Reagent 1	20 mL × 3	78234
				Reagent 2	6 mL × 3	78403
ChE	CicaFit ChE J	JSCC (pHBC)	7170	Reagent 1	60 mL × 4	78253
				Reagent 2	20 mL × 4	78254
LAP	CicaFit LAP	Leu-PNA	SPS	Reagent 1	20 mL × 3	77151
				Reagent 2	9 mL × 3	78404
<b>Blood Glucose</b>						
GLU	CicaLiquid GLU J	HK·G-6-PD	7170	Reagent 1	60 mL × 4	78224
				Reagent 2	20 mL × 4	77846
<b>Lipids</b>						
NEFA	CicaLiquid NEFA	Enzymatic	SPS	Reagent 1	20 mL × 3	77141
				Reagent 2	10 mL × 3	77142
<b>Protein</b>						
TP	CicaLiquid TP	Biuret	7170	Reagent 1	60 mL × 4	77831
				Reagent 2	38 mL × 4	78353
ALB	CicaLiquid ALB	BCG	7170	Reagent 1	60 mL × 4	77835
<b>Nitrogen</b>						
T-BIL	CicaLiquid T-BIL	Enzymatic	7170	Reagent 1	60 mL × 4	77801
				Reagent 2	19 mL × 4	78354
D-BIL	CicaLiquid D-BIL	Enzymatic	7170	Reagent 1	60 mL × 4	77803
				Reagent 2	19 mL × 4	78355
UA	CicaLiquid UA	Uricase	7170	Reagent 1	60 mL × 4	77476
				Reagent 2	20 mL × 4	77477
UN	CicaLiquid-N UN	Enzymatic	7170	Reagent 1	60 mL × 4	77679
				Reagent 2	20 mL × 4	77680
CRE	CicaLiquid-N CRE	Enzymatic	7170	Reagent 1	60 mL × 4	77663
				Reagent 2	20 mL × 4	77664
NH <sub>3</sub>	CicaLiquid NH <sub>3</sub>	Enzyme cycling	SPS	Reagent 1	20 mL × 3	77836
<b>Urine</b>						
m-TP	CicaLiquid m-TP	Pyrogallol Red	SP	Reagent 1	20 mL × 4	78207
NAG	CicaFit NAG	3MP-NAG	SPS	Reagent 1	20 mL × 3	78209
				Reagent 2	10 mL × 3	78210

Test Item	Product Name	Method	Type	Package	Code	
<b>Electrolites</b>						
IP	CicaLiquid IP	Enzymatic	7170	Reagent 1	60 mL × 4	77859
				Reagent 2	60 mL × 2	77860
Ca	CicaFit Ca	CPZIII	7170	Reagent 1	60 mL × 4	77227
				Reagent 2	33 mL × 4	78357
Mg	CicaLiquid Mg	Enzymatic	7170	Reagent 1	50 mL × 2	77193
				Reagent 2	20 mL × 2	77194
Fe	CicaLiquid Fe	Nitroso-PSAP	7170	Reagent 1	50 mL × 4	77626
				Reagent 2	20 mL × 4	77627
UIBC	CicaLiquid UIBC	Nitroso-PSAP	7170	Reagent 1	50 mL × 4	77630
				Reagent 2	20 mL × 4	77631

## Immunology

Test Item	Product Name	Method	Type	Package	Code	
<b>Immunology</b>						
CRP	Cias V CRP	Latex Turbidimetric Immunoassay	7170	Reagent 1	60 mL × 3	77950
				Reagent 2	60 mL × 3	77951
m-ALB	Cias ALB-M	Turbidimetric Immunoassay	—	Reagent 1	20 mL × 3	77427
				Reagent 2	10 mL × 2	
PG	Cias Latex PG I	Latex turbidimetry	S	Reagent 1	20 mL × 2	78227
	Cias Latex PG II	Latex turbidimetry		Reagent 2	8 mL × 1	
Cys-C	Cias Cys-C	Latex turbidimetry	7170	Reagent 1	20 mL × 2	78228
				Reagent 2	8 mL × 1	
MMP-3	Cias MMP-3	Latex turbidimetry	S	Reagent 1	15 mL × 2	78244
				Reagent 2	6 mL × 1	
H-FABP	Cias H-FABP	Latex turbidimetry	7170	Reagent 1	20 mL × 2	78269
				Reagent 2	10 mL × 1	
H-FABP	Cias H-FABP	Latex turbidimetry	7170	Reagent 1	20 mL × 2	78247
				Reagent 2	7 mL × 1	

## Calibrators & Controls

Test Item	Product Name	Package	Storage	Code
<b>Clinical Chemistry Calibrators</b>				
GLU	Glucose Calibrator	10 mL × 6	2~8°C	77094
TP, ALB	TP/ALB Calibrator	3 mL × 3	2~8°C	77811
T-BIL, D-BIL	BIL Calibrator	2 mL × 3	2~10°C	77696
UA	UA Calibrator	10 mL × 6	2~8°C	77100
UN	UN Calibrator	10 mL × 6	2~8°C	77096
CRE	CRE Calibrator	10 mL × 6	2~8°C	77097
NH <sub>3</sub>	NH <sub>3</sub> Calibrator	10 mL × 3 × 2 conc.	2~8°C	77850
IP	IP Calibrator	10 mL × 6	2~8°C	77099
Ca	Ca Calibrator	10 mL × 6	2~8°C	77098
Mg	Mg Calibrator	10 mL × 6	2~8°C	77101
Fe	Fe Calibrator	10 mL × 6	2~8°C	77285
NEFA	NEFA Calibrator	3 mL × 3	2~8°C	77172
m-TP	m-TP Calibrator	3 mL × 3	2~8°C	78208
NAG	NAG Calibrator (F)	1 mL × 3	Below -20°C	78311

Test Item	Product Name	Package	Storage	Code
<b>Multianalyte Calibrators</b>				
AST, ALT, $\gamma$ -GT, CK, LD, ALP, AMY, p-AMY, ChE	TraceCalib PLUS	1 mL $\times$ 3	-30 $\sim$ -20 $^{\circ}$ C	77939
CK-MB, CK	CK-MB Calib	1 mL $\times$ 3	Below -20 $^{\circ}$ C	78202
CRE, GLU, IP, Ca, Fe, Mg	MultiCalib 6M	5 mL $\times$ 3	2 $\sim$ 8 $^{\circ}$ C	77937
UA, UN, GLU	MultiCalib 3M	5 mL $\times$ 3	2 $\sim$ 8 $^{\circ}$ C	77938
AST, ALT, $\gamma$ -GT, CK, LD, ALP, AMY, p-AMY, ChE, LAP, GLU, CRE, UA, UN, T-BIL, D-BIL, TP, ALB, Fe, Ca, Mg, IP	CalibSera 24	1 mL $\times$ 3	Below -30 $^{\circ}$ C	78201
<b>Immunology Calibrators</b>				
CRP	V CRP Calibrator (A)	1 mL $\times$ 1 $\times$ 5 conc.	2 $\sim$ 8 $^{\circ}$ C	77953
m-ALB	Multi ALB-M Calibrator	2 mL $\times$ 1 $\times$ 5 conc.	2 $\sim$ 8 $^{\circ}$ C	77879
PG I, PG II	PG Calibrator (S)	3 mL $\times$ 1 (Blank) 1 mL $\times$ 1 $\times$ 5 conc.	2 $\sim$ 8 $^{\circ}$ C	78229
Cys-C	Cys-C Calibrator	1 mL $\times$ 1 $\times$ 6 conc.	2 $\sim$ 8 $^{\circ}$ C	78245
MMP-3	MMP-3 Calibrator (S)	1 mL $\times$ 1 $\times$ 6 conc.	2 $\sim$ 8 $^{\circ}$ C	78270
H-FABP	H-FABP Calibrator	1 mL $\times$ 1 $\times$ 6 conc.	2 $\sim$ 8 $^{\circ}$ C	78248
<b>Control Materials</b>				
AST, ALT, $\gamma$ -GT, CK, LD, ALP, AMY, p-AMY, ChE	TraceCheck PLUS	2 mL $\times$ 3 $\times$ 2 conc.	-30 $\sim$ -20 $^{\circ}$ C	77940
CK-MB	Qualitrol CK-MB (L)	1 mL $\times$ 3	Below -20 $^{\circ}$ C	78203
AST, ALT, $\gamma$ -GT, CK, LD, ALP, AMY, p-AMY, ChE, LAP, GLU, CRE, UA, UN, TP, ALB, T-BIL, D-BIL, Fe, Ca, Mg, IP, CRP	L-Qualitrol	3 mL $\times$ 3 $\times$ 2 conc.	Below -20 $^{\circ}$ C	78200
NH <sub>3</sub>	Ammonia Control	5 mL $\times$ 3 $\times$ 2 conc.	2 $\sim$ 8 $^{\circ}$ C	77878
NAG	Qualitrol NAG (F)	2 mL $\times$ 3	Below -20 $^{\circ}$ C	78312
PGI, PG II	Qualitrol PG (S)	2 mL $\times$ 2 $\times$ 2 conc.	2 $\sim$ 8 $^{\circ}$ C	78230
Cys-C	Qualitrol Cys-C	1 mL $\times$ 3 $\times$ 2 conc.	2 $\sim$ 8 $^{\circ}$ C	78246
MMP-3	Qualitrol MMP-3 (S)	1 mL $\times$ 3 $\times$ 2 conc.	2 $\sim$ 8 $^{\circ}$ C	78271
H-FABP	Qualitrol H-FABP	1 mL $\times$ 3 $\times$ 2 conc.	2 $\sim$ 8 $^{\circ}$ C	78249



**KANTO CHEMICAL CO., INC.**

**REAGENT DIVISION**

2-1, Nihonbashi Muromachi 2-chome, Chuo-ku, Tokyo, 103-0022, Japan

TEL : +81-3-6214-1092

E-mail : kanto-61@kanto.co.jp

Website: <https://www.kanto.co.jp>